

HEARTLIGHT® X3 - US ENDOSCOPIC ABLATION SYSTEM



Product Specifications

	HeartLight X3 Catheter	HeartLight Sheath	HeartLight Endoscope	Balloon Fill Media	HeartLight X3 Console
					
Product Code	18-5000	18-3356	18-1447	18-1342	18-4952
French Size	12F	12F (ID)	2F	N/A	N/A
Dimensions	340-360cm (110cm usable length)	90cm (75cm usable length)	300cm (145cm usable length)	N/A	61 x 76 x 160cm (W x H x D) 140 kg
Other	Recommended storage: 5°C - 30°C	Use only 0.035 or 0.038 extra stiff exchange length guidewires 180° unidirectional deflection	Refer to the Endoscope Cleaning and Sterilization Instructions (06-5123)	D2O and Sodium Diatrizoate (contrast media)	Power requirements: 120 V 50-60 Hz



CardioFocus, Inc.
500 Nickerson Road, Suite 500-200
Marlborough, Massachusetts 01752

T: (508) 658-7200
information@cardiofocus.com

